



Ingredient List

75ml water

10gr
black
cabbage

10gr
glycerol

20gr
gelatine

QR CODE

Preparation Steps

POOR

MIX


BOIL 5
MIN

POUR THE
MIXTURE IN
MOLD

Material Assessments

hard	<input type="checkbox"/>	soft	<input type="checkbox"/>	opaque	<input type="checkbox"/>	transparent	<input type="checkbox"/>
smooth	<input type="checkbox"/>	rough	<input type="checkbox"/>	tough	<input type="checkbox"/>	ductile	<input type="checkbox"/>
matte	<input type="checkbox"/>	glassy	<input type="checkbox"/>	light	<input type="checkbox"/>	heavy	<input type="checkbox"/>
elastic	<input type="checkbox"/>	not elastic	<input type="checkbox"/>	fibred	<input type="checkbox"/>	not-fibred	<input type="checkbox"/>
strong	<input type="checkbox"/>	weak	<input type="checkbox"/>				

Forming Technique



References

WWW.MATERIOM.ORG



BIOPLASTIC LIBRARY

Computation Based Basic Design @ BILGI
SPRING 2022 - 23

ÜZEYİR DURMUŞ/AHMET HASAN KAYA