



Ingredient List

- 1ml glycerol
- 50ml water
- 3gr agar-agar-
- 10gr coffee




Preperation Steps

- POOR
- MIX
- BOIL 5 MIN
- POUR THE MIXTURE IN MOLD

Material Assessments

hard	<input type="checkbox"/>	soft	<input type="checkbox"/>	opaque	<input type="checkbox"/>	transparent	<input type="checkbox"/>
smooth	<input type="checkbox"/>	rough	<input type="checkbox"/>	tough	<input type="checkbox"/>	ductile	<input type="checkbox"/>
matte	<input type="checkbox"/>	glassy	<input type="checkbox"/>	light	<input type="checkbox"/>	heavy	<input type="checkbox"/>
elastic	<input type="checkbox"/>	not elastic	<input type="checkbox"/>	fibred	<input type="checkbox"/>	not-fibred	<input type="checkbox"/>
strong	<input type="checkbox"/>	weak	<input type="checkbox"/>				

Forming Technique



References

WWW.MATERIOM.ORG

