


VINEGAR/ STRACH
BIOMETARIAL



5ML GLYCEROL



5 ML VINEGAR



15 GR STRACH



65 ML WATER


Preperation Steps



Material Assessments

hard	• • • •	soft		opaque	• • • •	transparent	
smooth	• • •	rough		tough	• • • •	ductile	
matte	• • •	glossy		light	•	heavy	
elastic	• • •	not elastic		fibred	•	not-fibred	
strong	• • •	weak					

Forming Technique



References

<https://yjad.yildiz.edu.tr/storage/upload/pdfs/1656575220-tr.pdf>